

TITLE: Lab 14- Respiratory Physiology

PURPOSE:

The purpose of this lab is to use a spirometer. This will determine the lung capacity of a tidal volume, vital capacity, inspiratory capacity, inspiratory reserve volume, expiratory capacity, and expiratory reserve volume. The recording will be obtained from spirometer may be analyzed to determine the relative pulmonary condition of humans. The timed vital capacity (TVC) or forced expiratory volume (FEVT) will also be calculated. The impedance pneumography and the measurement of ventilation rates will be recording the changing impedance of an expanding and contracting thorax.

PROCEDURE:


14-B


1. The Morgan ComPAS computer program has already calculated and factored in the BTPS (Body Temperature Pressure Saturation) correction factor.
2. Fully insert the Pneumotrac filter/mouthpiece you purchased at the bookstore. If you have difficulty keeping air from leaking through your nose, you may need to wear a nose clip, as air leakage will result in inaccurate results.
3. Be sure the correct student information is loaded up before you start the FVC test.
4. After starting the FVC test, follow the verbal instructions of your instructor: begin with your mouth off the mouthpiece so the pneumotach can equilibrate; after getting a good seal with your mouth, start with tidal breathing; when you are ready, take in the deepest breath possible, then forcefully blow it out as fast as you can and keep squeezing until instructed to stop. The instructor will print out your "FVC Volume Time Curve" (part of your 14-B results), and it should look similar to Figure 14-2
5. To calculate the vital capacity for the FVC test (also called the forced expiratory volume), measure the height of the highest peak of the curve in mm and multiply that length in mm by 66.67 ml/mm (our FVC conversion factor). Then round off ml to whole numbers. (NOTE: this is similar to the 14-A SVC calculations, but with a different conversion factor.) Just like in 14-A, use the gridlines to double check that your figures are in the ballpark (e.g., if you calculated the vital capacity in Fig. 14-2 on p. 94 to be 3635 ml, you must be off because you can tell from just looking at the gridlines that it is much closer to 4500 ml than 3635 ml). Can you see this in Fig. 14-2?
6. Go to the "1 second" vertical line in your FVC graph and measure the height where the curved line crosses the 1 second vertical line in the same way as you did for the FVC in step 5. This is your FEV1 volume.
7. Divide the volume you calculated for FEV1 by the volume you calculated for the vital capacity in step 6, and then multiply by 100 to determine the percentage of the vital capacity exhaled at one second.
8. Go to the "3 second" vertical line in your FVC graph and measure the height where the curved line crosses the 3 second vertical line in the same way as you did in steps 5 and 6. This is your FEV3 volume.
9. Divide the volume you calculated for FEV3 by the volume you calculated for the vital capacity in step 6, and then multiply by 100 to determine the percentage of the vital capacity exhaled at three seconds.
10. Compare these values to the predicted values and explain possible causes for any differences.

14-C

1. Open the grey plastic box on your lab desk that says “BASELINE LungCapacitySpirometer” on the lid. Inside the lid of the box is a white paper that has specific instructions, please read the whole inside page with “how to use.”
2. Insert the clear plastic mouthpiece on the “Windmill-Type” spirometer and make sure the measurement indicator is at the zero position before beginning.
3. Make sure you only exhale into the spirometer, DO NOT inhale from it.
4. After exhaling, record the measurement from the spirometer. Be sure to place your used plastic mouthpiece in the correct tube after use (the tube is labeled).
5. Calculate your predicted vital capacity from the nomograms available in lab. Using a straightedge, make a line matching your height and age to the vital capacity prediction. Note that the VC is in liters whereas other measurements have been taken in milliliters.
6. Compare the values obtained from the portable spirometer, the predicted values from the nomograms, and the value obtained from the Koko spirometer, if available. How can you account for any differences? (NOTE: your predicted VC from the nomogram, and a comparison to the measured VC in 14-A should be included in your discussion of 14-A)

RESULTS:





DISSCTION:

There's multiple lung volumes and capacities are measures of the amount of an air in the respiratory system. Some of them are tidal volume (TV) the average is 500mL, inspiratory reserve volume (IRV) the average is around 3,000mL. The maximum volume of air that can be inhaled after normal tidal volume inhalation. There's expiratory reserve volume (ERV) the average is approximately 1,200mL. The maximum volume of air that can be exhaled after a normal tidal volume exhalation. There's also residual volume (RV) the average is 1,200mL. The volume of air remaining in the lungs after a maximal exhalation. There's vital capacity (VC) the average is 4,500mL. It's the maximum volume of air that can exhaled after maximal inhalation it's the sum of tidal volume, inspiratory reserve volume and expiratory reserve volume. There's also total lung capacity (TLC) the average is approximately 5,700mL. The maximum volume of the air the lungs can hold its sum of the vital capacity and residual volume. Then, there's functional residual capacity (FRC) the average is around 2,400mL. The volume of air remaining in the lungs after a normal tidal volume. Finally, this some functions of the lung capacity and there average values providing general reference.

CONCLUSION:

In conclusion, the impedance pneumography is a non-invasive method used to measure respiratory parameters. The respiratory rate is calculated by measuring the frequency of impedance changes associated with each breath. The tidal volume is indirectly estimated by assessing the amplitude of impedance changes. There's also minute ventilation is calculated by multiplying respiratory rate by tidal volume. The inspiratory time and expiratory time these times can be derived from the shape and duration of impedance waveform. The derived parameters can be system may be provided derived parameter such as dynamic lung compliance and other indices related to respiratory mechanic. These are important to know impedance pneumography is a valuable tool for assisting respiratory function